




# QuickCure



## Marcy QuickCure Sensor Gel Refill


Marcy's QuickCure gel saves you both time and money. It cures real quick—half the time compared to other gels. And the QC12, with three times the gel (and extra nozzles), gives you multiple uses with a single dispenser. This cuts waste and cost. Quick, convenient, and cost-saving. Marcy QuickCure.

- **Save Time**—cures in under 4 minutes at 70° (less than half the time of leading gels)
- **Cut Waste**—use just what you need with the multi-use QC12
- **Save Money**—lower price, plus quick cure and no-waste means big savings
- **Clean & easy removability when fully cured**—great for repeat customer/warranty work

Part No.	Description
<b>QC4</b>	4ml (2 x 2ml) syringe with mixing nozzle
<b>QC12</b>	12ml (2 x 6ml) syringe with 3 mixing nozzles
<b>QCN</b>	Mixing nozzle replacements (x10)

## Marcy QuickCure vs. other gel

### SENSOR GEL CURING TIMES AT DIFFERENT TEMPERATURES



\*To replicate the most demanding conditions, Marcy QuickCure gel was dispensed only after it reached ambient temperature. Quicker cure times can be achieved, even in sub-freezing temperatures, by keeping the gel at room temperatures prior to dispensing.

Available at your local distributor—contact Marcy for details

1-614-471-5200 | [marcyadhesives.com](http://marcyadhesives.com)